


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
o-Cresol - ReagentPlus®, ≥99%

Product Number: **C85700**
CAS Number: 95-48-7
MDL: MFCD00002226
Formula: C₇H₈O
Formula Weight: 108.14 g/mol



TEST	Specification
Appearance (Color) Colorless to Light Yellow (Liquid Phase) Colorless to White to Off-White (Solid Phase)	Conforms to Requirements
Appearance (Form) Liquid or Granules, Crystals, Crystalline Solid, or Solid	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 99.0 %
Solubility (Turbidity) Solubility Concentration: 100 mg/mL in 95% Ethanol	Clear
Solubility (Color) Solubility Concentration: 100 mg/mL in 95% Ethanol	Colorless to Faint Yellow

Specification: PRD.2.ZQ5.10000003539

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.