

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Phosphatase, Acid from potato - lyophilized powder, ≥ 3.0 units/mg solid

Product Number:

P1146

CAS Number:

9001-77-8

MDL:

MFCD00131847

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

units/mg solid

 ≥ 3.0 One unit will hydrolyze 1.0 micromole of
p-nitrophenyl phosphate per minute at pH 4.8 at 37 deg C.

Recommended Retest Period

2 Years

Specification: PRD.2.ZQ5.10000030890

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.