


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Hydrocinnamic acid - 99%

Product Number: 135232
CAS Number: 501-52-0
MDL: MFCD00002771
Formula: C₉H₁₀O₂
Formula Weight: 150.17 g/mol



TEST

Specification

Appearance (Color)	White to Off White
Appearance (Form)	Powder or Crystal or Chunk(s)
Infrared spectrum	Conforms to Structure
Purity (Titration by NaOH)	98.5 - 101.5 %
Purity (GC)	≥ 98.5 %
Solubility (Turbidity)	Clear
50mg/ml, Ethanol	
Solubility (Color)	Colorless to Faint Yellow

Specification: PRD.0.ZQ5.10000035130

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.