

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Silver – Nanowires, diam. × L 50 nm (±10 nm) × 40 μm (±5 μm), 5 mg/mL (in IPA)

Product Number:

807508

Storage Temperature:

2 - 8 °C

TEST	Specification
Appearance (Color)	Grey
Appearance (Form)	Suspension
ICP Major Analysis Confirms Silver Component	Confirmed
Concentration 5 mg/mL	Confirmed
Diameter 50+/-10 nm	Confirmed
Length 40+/-5 μm	Confirmed
Isopropyl Alcohol	Confirmed
Recommended Retest Period 1 Year	-----

Specification: PRD.0.ZQ5.10000096590

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.