


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Palau'Chlor® - 95%

Product Number: 792454

Formula: C₅H₈ClN₃O₄
Formula Weight: 209.59 g/mol
Storage Temperature: 2 - 8 °C**TEST****Specification**

Appearance (Color)	White to Beige
Appearance (Form)	Powder or Crystals
Proton NMR spectra	Conforms to Structure
Titration with AgNO ₃ After Oxygen Combustion	94.5 - 105.5 %
Registered Trademark Palau'Chlor is a registered trademark of Sigma-Aldrich Co. LLC	Confirmed

Specification: PRD.1.ZQ5.10000047939

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.