

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

1,1,2,2-Tetrakis(p-hydroxyphenyl)ethane - powder

Product Number:

H1259

CAS Number:

7727-33-5

MDL:


MFCD00002371

Formula:

C<sub>26</sub>H<sub>22</sub>O<sub>4</sub>

Formula Weight:

398.45 g/mol



### TEST

### Specification

Appearance (Color)

White to Off-White

Appearance (Form)

Powder

Solubility (Color)

Colorless to Light Yellow

Solubility (Turbidity)

Clear

50 mg/mL, MeOH

Proton NMR spectrum

Conforms to Structure

% Solvent

≤ 10

by GC

Specification: PRD.2.ZQ5.10000012543

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.