

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Oleic Acid-Water Soluble - powder, BioReagent, suitable for cell culture

Product Number: **O1257**

MDL: MFCD00064242



Storage Temperature: -20 °C

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear to Slightly Hazy
250 mg + 5 mL H ₂ O	
Cell Culture Test	Pass
mg/g	≥ 28
by HPLC	
mg oleic acid/g solid	

Specification: PRD.0.ZQ5.10000039603

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.