


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ethyl diazoacetate solution - 15% in toluene

Product Number: 752150



Formula: C₄H₆N₂O₂
Formula Weight: 114.10 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	Yellow to Dark Yellow and Orange to Dark Orange and Brown to Dark Brown
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Proton NMR spectrum	Conforms to Structure
Weight wt% of Ethyl Diazoacetate	14.0 - 16.5

Specification: PRD.1.ZQ5.10000015996

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.