


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Thymoquinone - ≥98%

Product Number: 274666
CAS Number: 490-91-5
MDL: MFCD00001602
Formula: C₁₀H₁₂O₂
Formula Weight: 164.20 g/mol



TEST

Specification

Appearance (Color)

White to Dark Yellow and Faint Orange to Orange and Faint Brown to Brown and Orange-Brown and Brown-Orange

Appearance (Form)

Powder, Crystals or Crystalline Powder and/or Chunks

Conforms to Requirements

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 98.0 %

Specification: PRD.1.ZQ5.10000049072

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.