

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Rabbit IgG (whole molecule)-Alkaline Phosphatase antibody produced in goat - IgG fraction of antiserum

Product Number:

A0418

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Almost Colorless to Light Yellow
IgG Concentration (OD 280nm)	0.8 - 2.5 mg/ml
ELISA - Titer	1:20,000-1:30,000
WB-Cell Line/Tissue Extract	Hela Cells
WB Titer	1:100,000-1:200,000
IHC - Tissue	Human Tonsil
IHC - Fixation	FFPE-Formalin Fixed Paraffin Embedded
IHC-Titer	1:40-1:80
Recommended Retest Period	-----
2 Years	

Specification: PRE.0.ZQ5.10000018119

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.