


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Heptanes, mixture of isomers - contains $\geq 25\%$ n-heptane, ReagentPlus*

Product Number: 730491



Formula: C₇H₁₆
Formula Weight: 100.20 g/mol

TEST	Specification
Color	\leq 10 APHA
GC Purity (sum of isomers)	\geq 99.0 %A
Assay	\geq 25 %
n-Heptane by GC	
Water (by Karl Fischer)	\leq 0.01 %
Boiling Range	Pass
92.0°C - 100.0°C	
Specific Gravity	0.690 - 0.720
@ 15.6 C	
Residue After Evaporation	\leq 0.001 %
Bromine Index	\leq 25
mg Br per 100 grams	

Specification: PRD.3.ZQ5.10000010830

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.