

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tellurium - powder, -30 mesh, 99.997% trace metals basis

Product Number: 264865
CAS Number: 13494-80-9
MDL: MFCD00134062
Formula: Te
Formula Weight: 127.60 g/mol

Te

TEST	Specification
Appearance (Color) Grey to Dark Grey	Conforms to Requirements
Appearance (Form) Powder or Granular Powder	Conforms to Requirements
Particle Size -30 Mesh	Conforms
X-Ray Diffraction	Conforms to Structure
Purity 99.997% Based On Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	≤ 35.0 ppm

Specification: PRD.1.ZQ5.10000007996

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.