


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Methylpyrazole-5-carboxylic acid - 97%

Product Number: 644927
CAS Number: 402-61-9
MDL: MFCD00462235
Formula: C₅H₆N₂O₂
Formula Weight: 126.11 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Powder, Crystals, Crystalline Powder and/or Chunks
Infrared spectrum
Purity (Titration by NaOH)
Purity (GC)

White to Tan
Conforms to Requirements

Conforms to Structure
96.5 - 103.5 %
≥ 96.5 %

Specification: PRD.0.ZQ5.10000012974

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.