

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Estrogen-Related Receptor α antibody produced in rabbit - affinity isolated antibody, buffered aqueous solutionProduct Number: **E0406**

MDL: MFCD06798438

Storage Temperature: -70 °C

TEST	Specification
Formulation	Conforms
Solution in phosphate buffered saline, and 0.1% sodium azide	
Concentration	Conforms
1 mg/mL	
Working Dilution	Conforms
by Immunohistology	
Recommended working dilution: 18 - 24 ug/mL using human tissue. Optimal dilution to be determined by the researcher.	

Specification: PRD.1.ZQ5.10000041536

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.