


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
5-Bromo-2-fluorophenylboronic acid

Product Number: 593621
CAS Number: 112204-57-6
MDL: MFCD05664307
Formula: C₆H₅BBrFO₂
Formula Weight: 218.82 g/mol



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Proton NMR spectrum	Conforms to Structure
Titration with NaOH	≥ 95.0 %

Specification: PRD.0.ZQ5.10000037577

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.