

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Fluoroindole - 97%

Product Number: 457396
CAS Number: 387-43-9
MDL: MFCD00055992
Formula: C₈H₆FN
Formula Weight: 135.14 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
White to Pink to Brown	
Appearance (Form)	Liquid or Solid
Infrared spectrum	Conforms to Structure
Carbon	68.6 - 73.6 %
Nitrogen	9.9 - 10.7 %
Purity (TLC)	≥ 97 %
Solubility (Turbidity)	Clear
50mg/ml, CHCl ₃	
Solubility (Color)	Conforms to Requirements
Colorless to Pink to Green to Brown	

Specification: PRD.0.ZQ5.10000032295

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.