

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2,3,4,6-Tetrafluorophenylboronic acid - ≥95%

Product Number:

680621

CAS Number:

511295-00-4

MDL:


MFCD04039309

Formula:

C₆H₃BF₄O₂

Formula Weight:

193.89 g/mol



TEST

Specification

Appearance (Color)

White to Tan

Appearance (Form)

Conforms to Requirements

Solid, Powder, Crystals, Crystalline Powder, and/or Chunks

Proton NMR Spectrum

Conforms to Structure

Titration with NaOH

≥ 95.0 %

Specification: PRD.0.ZQ5.10000037997

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.