

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Hoveyda-Grubbs Catalyst™ 2nd Generation - 97%

Product Number: 569755
CAS Number: 301224-40-8
MDL: MFCD03701614
Formula: C₃₁H₃₈Cl₂N₂ORu
Formula Weight: 626.62 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Conforms to Requirements
Green	
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Carbon	57.3 - 61.5 %
Nitrogen	4.2 - 4.7 %

Specification: PRD.1.ZQ5.1000005404

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.