

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Cellulase from *Trichoderma viride* - crude powder, 3-10 units/mg solid

Product Number:

C9422

CAS Number:

9012-54-8

MDL:

MFCD00081510

Storage Temperature:

2 - 8 °C

TEST

Specification

units/mg solid

3 - 10

Enzymatic Activity: One unit will liberate 1.0 micromole of glucose from cellulose per hour at pH 5.0 at 37 deg C (2 hour assay)

Specification: PRD.0.ZQ5.10000041432

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.