

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Platinum on alumina - extent of labeling: 5 wt. % loading, -325 mesh, powder

Product Number:

311324

Pt

MDL:

MFCD00011179

| TEST | Specification |
|------|---------------|
|------|---------------|

Appearance (Color)

Conforms to Requirements

Dark Grey to Brown

Appearance (Form)

Powder

Loss on Drying

≤ 5.0 %

Size

Conforms

-325 Mesh

ASSAY

4.0 - 6.0 %

Platinum by ICP

Specification: PRD.0.ZQ5.10000029487

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.