

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Tungsten hexacarbonyl - 99.99% trace metals basis (excluding Mo), purified by sublimation

**Product Number:** 472956  
**CAS Number:** 14040-11-0  
**MDL:** MFCD00011462  
**Formula:** C<sub>6</sub>O<sub>6</sub>W  
**Formula Weight:** 351.90 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Carbon Content	20.17 - 20.79 %
ICP Major Analysis	Confirmed
Confirms Tungsten Component	
Purity	Meets Requirements
99.99% Based on Trace Metals Analysis Excluding Mo	
Trace Metal Analysis	≤ 150.0 ppm

Specification: PRD.1.ZQ5.10000039374

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.