

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Silver(I) fluoride - $\geq 99.9\%$ trace metals basis

Product Number: 226858
CAS Number: 7775-41-9
MDL: MFCD00003410
Formula: AgF
Formula Weight: 126.87 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Yellow to Very Dark Yellow and Orange to Very Dark Orange and Brown to Very Dark Brown and Yellow-Orange and Yellow-Brown and Orange-Brown
Appearance (Form)	Powder
Titration with KSCN % Silver	83.7 - 86.3 %
ICP Major Analysis Confirms Silver Component	Confirmed
Purity > = 99.9% Based On Trace Metals Analysis	Conforms
Trace Metal Analysis	≤ 1000.0 ppm

Specification: PRD.1.ZQ5.10000039097

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.