

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Formaldehyde Dehydrogenase from Pseudomonas sp. - lyophilized powder, 1.0-6.0 units/mg solid

Product Number: F1879
CAS Number: 9028-84-6
MDL: MFCD00131100

Storage Temperature: -20 °C

TEST**Specification**

Appearance (Color)	White to Off White
Appearance (Form)	Powder or Crystals
units/mg solid	1.0 - 6.0
One unit will oxidize 1.0 micromole of formaldehyde to formic acid per minute at pH 7.5 at 37 deg C	

Specification: PRD.2.ZQ5.10000016859

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.