

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-CCL17/TARC antibody produced in goat - affinity isolated antibody, lyophilized powder

Product Number:

C1497

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Form) from a 0.2 um filtered solution in PBS with 5% trehalose	Lyophilized Powder
Endotoxin Assay	Less than 0.1 Eu per 1 ug
Western Blot 0.1 - 0.2 µg/ml	Conforms
Immunohistochemistry 5 - 15 µg/ml	Conforms

Specification: PRD.1.ZQ5.10000012144

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.