


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(-)- α -Pinene - $\geq 97\%$, FCC

Product Number: W290203
CAS Number: 7785-26-4
MDL: MFCD00064145
Formula: C₁₀H₁₆
Formula Weight: 136.23 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.464 - 1.468
Infrared spectrum	Conforms to Structure
Purity (GC)	$\geq 97.0 \%$
Specific Gravity	0.855 - 0.860
Solubility (Turbidity)	Clear
Solubility (Color)	Colorless
1ml/3ml,95% Ethanol	
Angular Rotation	-50.0 - -20.0 °
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
5 Years	

Specification: PRD.1.ZQ5.10000034555

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.