


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Perfluorohexyl iodide - 99%

Product Number: 262145
CAS Number: 355-43-1
MDL: MFCD00001063
Formula: C₆F₁₃I
Formula Weight: 445.95 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

Colorless to Pink to Orange or Purple

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

GC (area %)

≥ 98.5 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.