


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
5-Methoxyindole - 99%

Product Number: M14900
CAS Number: 1006-94-6
MDL: MFCD00005674
Formula: C₉H₉NO
Formula Weight: 147.17 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Grey and Faint Beige to Beige and Faint Brown to Brown and Grey-Brown
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Carbon	72.3 - 74.6 %
Nitrogen	9.1 - 9.9 %
Purity (GC)	≥ 98.5 %

Specification: PRD.0.ZQ5.10000011763

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.