

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Zinc acrylate - 98%

Product Number: 409464
CAS Number: 14643-87-9
MDL: MFCD00067208
Formula: C₆H₆O₄Zn
Formula Weight: 207.50 g/mol

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Complexometric EDTA (% Zn)	30.7 - 32.3 %

Specification: PRD.0.ZQ5.10000029454

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.