

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Methyl isocyanoacetate – technical grade, 95%

Product Number:**238880**

CAS Number:

39687-95-1

MDL:

MFCD00000006

Formula:


C₄H₅NO₂

Formula Weight:

99.09 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

Yellow to Dark Yellow and Orange to Dark
Orange and Brown to Dark Brown and Red to Dark Red and
Yellow-Orange and Yellow-Brown and Orange-Brown and Red-Brown

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 94.0 %

Specification: PRD.0.ZQ5.10000019586

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.