

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Platinum(IV) chloride - 96%

Product Number: 206113
CAS Number: 13454-96-1
MDL: MFCD00011182
Formula: Cl₄Pt
Formula Weight: 336.89 g/mol

PtCl₄

TEST	Specification
------	---------------

Appearance (Color)	Faint Brown to Very Dark Brown
Appearance (Form)	Powder or Chunks
Platinum (Pt)	55.3 - 60.5 %
ICP Major Analysis	Confirmed
Confirms Pt Component	

Specification: PRD.0.ZQ5.10000035788

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.