

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Amino-6-fluoropyridine - 98%

Product Number: 725064
CAS Number: 1597-32-6

Formula: C₅H₅FN₂
Formula Weight: 112.11 g/mol



TEST

Specification

Appearance (Color)

White to Grey and Faint Beige to Beige
and Faint Brown to Brown and Grey-Brown

Appearance (Form)

Powder or Crystals

Infrared spectrum

Conforms to Structure

Purity (Titration by HClO₄)

97.5 - 102.5 %

Purity (GC)

≥ 97.5 %

Recommended Retest Period

2 Years

Specification: PRD.1.ZQ5.10000008161

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.