


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Perylene - sublimed grade, ≥99.5%

Product Number: 394475
CAS Number: 198-55-0
MDL: MFCD00004142
Formula: C₂₀H₁₂
Formula Weight: 252.31 g/mol



TEST

Specification

Appearance (Color)	Orange
Appearance (Form)	Powder or Crystal or Chunk(s)
Infrared spectrum	Conforms to Structure
Melting point by DSC	
Thermogravimetric Analysis	
Report Temperature at 0.5% Weight Loss	
Purity (HPLC)	≥ 99.50 %

Specification: PRD.1.ZQ5.10000030041

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.