

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Platinum – nanopowder, 200 nm particle size (SEM), 99.9% (metals basis)

Product Number: 771937 Pt
CAS Number: 7440-06-4
Formula: Pt
Formula Weight: 195.08 g/mol

TEST	Specification
Appearance (Color)	Grey
Appearance (Form)	Powder
Particle Size	Conforms
200 nm (SEM)	
ICP Major Analysis	Confirmed
Confirms Platinum Component	
Trace Metal Analysis	≤ 1500.0 ppm
Purity	Meets Requirements
99.9% Based On Trace Metals Analysis	

Specification: PRD.0.ZQ5.10000038968

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.