


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
XtalFluor-M<sup>®</sup>Product Number: **719447**  
CAS Number: 63517-33-9Formula: C<sub>4</sub>H<sub>8</sub>BF<sub>6</sub>NOS  
Formula Weight: 242.98 g/mol  
Storage Temperature: -20 °C

### TEST

### Specification

Appearance (Color)  
Appearance (Form)  
Infrared spectrum  
Proton NMR spectrum  
Fluorine NMR SpectraWhite to Beige and Faint Brown to Brown  
Powder or Solid  
Conforms to Structure  
Conforms to Structure  
Conforms to Structure

Specification: PRD.1.ZQ5.10000038862

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.