

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Malonic acid - ReagentPlus*, 99%

Product Number: M1296 CAS Number: 141-82-2

MDL: MFCD00002707

Formula: C3H4O4 Formula Weight: 104.06 g/mol HO OF

TEST

Appearance (Color)

Appearance (Form) Infrared spectrum Purity (Titration by NAOH) GC (area %)

Solubility (Turbidity)

c = 100 mg/ml; 1N NaOH

Solubility (Color)

Specification

White

Powder or Crystals or Granules or Chunks Conforms to Structure 98.5 - 101.5 % > 98.5 % Clear to Slightly Hazy

Colorless to Faint Yellow

Specification: PRD.0.ZQ5.10000028318

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.