

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
SODIUM NITRITE, ANHYDROUS, REDI-DRI (TM&

Product Number: 901903
CAS Number: 7632-00-0

Formula: NaNO₂
Formula Weight: 69.00 g/mol

TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Turbidity)	Clear
Solubility (Color)	Colorless to Faint Yellow
Purity	≥ 99.0 %
Recommended Retest Period	-----
2 Years	

Specification: PRD.0.ZQ5.10000088758

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.