

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Alginic acid sodium salt from brown algae - BioReagent, plant cell culture tested, low viscosity, powder

Product Number: A0682
CAS Number: 9005-38-3
MDL: MFCD00081310

TEST	Specification
Appearance (Color)	White to Beige to Brown
Appearance (Form)	Powder
Solubility (Color)	Faint Yellow
Solubility (Turbidity)	Slightly Hazy to Hazy
10 mg/mL, H ₂ O	
Brookfield Viscosity	4 - 12 cps
(1% in H ₂ O at 25 deg C)	
Plant Cell Culture Test	Pass

Specification: PRD.3.ZQ5.10000010662

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.