

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

3-Aminopropionitrile fumarate salt - metabolite

Product Number:

A3134

CAS Number:

2079-89-2

MDL:

MFCD00078321

Formula:

C₃H₆N₂ · 0.5C₄H₄O₄

Formula Weight:

128.13 g/mol



TEST

Specification

Appearance (Color)

White to Faint Yellow

Appearance (Form)

Powder

Solubility (Color)

Colorless to Very Faint Yellow

Solubility (Turbidity)

Clear

50 mg/mL, H₂O

1H NMR Spectrum

Conforms to Structure

Carbon

45.9 - 47.8 %

Nitrogen

21.4 - 22.3 %

Purity (TLC)

≥ 98 %

Specification: PRD.3.ZQ5.10000038171

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

