

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Manganese(II) fluoride - 98%

Product Number: 339288
CAS Number: 7782-64-1
MDL: MFCD00016222
Formula: F₂Mn
Formula Weight: 92.93 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Conforms to Requirements
Faint Pink to Purple	
Appearance (Form)	Powder
ICP Major Analysis	Confirmed
Confirms Manganese Component	
Complexometric EDTA	57.6 - 60.6 %
% Mn	

Specification: PRD.0.ZQ5.10000032414

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.