

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com Email USA: techsery@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

ProteoMass™ ACTH Fragment 18-39 MALDI-MS Standard - vial of 10 nmol, monoisotopic mol wt 2,464.1989 Da

Product Number:

A8346 CAS Number: 53917-42-3 MDL: MFCD00076159 C112H165N27O36 Formula: Formula Weight: 2,465.7 g/mol

Manning and a

TEST	Specification

Protein (nanomoles/vial) 8.0 - 18.0

by Micro BCA µg/vial 19.7 - 44.4

by Micro BCA

Suitability Suitable

Suitable as a standard for MALDI-Mass Spectrometry

Specification: PRD.1.ZQ5.10000063876

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.