

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Reagent Alcohol - reagent grade

Product Number: 362808
CAS Number: 64-17-5
MDL: MFCD00003568
Formula: C₂H₆O
Formula Weight: 46.07 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared Spectrum	Conforms to Structure
Purity (GC)	≥ 89.0 %
Ethanol	
GC Identified Peak	3.0 - 5.0 %
Methanol	
GC Identified Peak	4.0 - 6.0 %
2-Propanol	
Water (by Karl Fischer)	≤ 0.10 %
Residue on Evaporation	≤ 0.0005 %

Specification: PRD.3.ZQ5.10000007651

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.