

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Monoclonal Anti-Cytokeratin Peptide 19 antibody produced in mouse – clone A53-B/A2, tissue culture supernatant

Product Number: C6930

Storage Temperature: -20 °C

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless to Faint Yellow
Protein by Biuret	5 - 20 mg/ml
Mouse IgG2a by RID	0.05 - 0.20 mg/ml
IHC - Tissue	Human Placenta
IHC - Fixation	FFPE-Formalin Fixed Paraffin Embedded
IHC - Titer	1:50-1:100
Recommended Retest Period	-----
5 Years	

Specification: PRE.0.ZQ5.10000016221

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.