

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

2-tert-Butyl-1,4-benzoquinone 98%

Product Number:

429074

CAS Number:

3602-55-9

MDL:


MFCD00666928

Formula:

C₁₀H₁₂O₂

Formula Weight:

164.2 g/mol

**TEST****Specification**

Appearance (Color) Yellow to Very Dark Yellow	Conforms to Requirements
Appearance (Form) Powder or Crystals	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Purity (GC)	≥97.5 %

Remarks:

Specification Date : 06/21/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.