

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Platinum – wire, diam. 0.076 mm, ≥99.99% trace metals basis

**Product Number:** 357375 Pt  
**CAS Number:** 7440-06-4  
**MDL:** MFCD00011179  
**Formula:** Pt  
**Formula Weight:** 195.08 g/mol

TEST	Specification
Appearance (Color) Silver	Conforms
Appearance (Form) Wire	Conforms
Thickness 0.076 mm	Conforms
ICP Major Analysis Confirms Pt Component	Confirmed
Trace Metal Analysis	≤ 100.0 ppm
Purity ≥99.99% Based On Trace Metals Analysis	Meets Requirements
Use Test Suitable for Thermocouple Use	Confirmed

Specification: PRD.1.ZQ5.10000038213

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.