

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium hydroxide – reagent grade, 90%, flakes

Product Number:	484016	KOH
CAS Number:	1310-58-3	
MDL:	MFCD00003553	
Formula:	HKO	
Formula Weight:	56.11 g/mol	

TEST	Specification
------	---------------

Appearance (Color)

White

Appearance (Form)

Conforms to Requirements

Flakes or Chunks

Titration by HCl

89.5 - 110.5 %

X-Ray Diffraction

Conforms to Structure

Specification: PRD.1.ZQ5.10000008616

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.