

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

C Reactive Protein from human fluids - buffered aqueous solution

Product Number: C4063
CAS Number: 9007-41-4
MDL: MFCD00130690

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Solution
Appearance (Turbidity)	Clear to Hazy
mg protein/ml (Lowry w/TCA)	≥ 0.5
mg protein/ml (Lowry w/TCA)	
Identity	Human Origin
By OIMM	
Purity by SDS Electrophoresis	≥ 90 %
Tested For Infectious Agents	Tested

Specification: PRD.1.ZQ5.10000013233

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.