


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Manganese(II) nitrate hydrate - 98%

Product Number: 288640
CAS Number: 15710-66-4
MDL: MFCD00149788
Formula: MnN2O6 · xH2O
Formula Weight: 178.95 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Pink
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Solid or Liquid	
ICP Major Analysis	Confirmed
Confirms Manganese Component	
Compleximetric EDTA	18.0 - 22.0 %
% Manganese	

Specification: PRD.0.ZQ5.10000028604

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.