

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Mineral oil - suitable for preparation of Nujol mulls for infrared spectroscopy, light oil

Product Number:

M3516

CAS Number:

8042-47-5

MDL:

MFCD00131611

TEST	Specification
------	---------------

Appearance (Turbidity)

Clear

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Brookfield Viscosity

≤ 30.0 cps

c = Neat at 25 Degrees Celsius

Specific Gravity

0.82 - 0.88

Suitability

Suitable

Suitable for Preparation of Nujol Mulls

Specification: PRD.0.ZQ5.10000027788

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.