

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Human IgG (Fc specific)-Cy3 antibody produced in goat - affinity isolated antibody, buffered aqueous solution

Product Number:

C2571

MDL:

MFCD00212909

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Pink
Appearance (Form)	Liquid
Protein Content	≥ 1.0 mg/ml
Conjugate concentration prior to the addition of BSA	
F/P Molar Ratio	3 - 9
Conjugate concentration prior to the addition of BSA	
Direct Histology Titer	Minimum Dilution 1:30
Tissue: Paraffin-embedded human tonsil	

Specification: PRD.0.ZQ5.10000017552

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.