

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Anti-MAP Kinase (ERK-1, ERK-2) antibody produced in rabbit – whole antiserum

Product Number: **M5670**

Storage Temperature: -20 °C

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Very Light Brown to Very Dark Brown and Dark Brown-Red to Very Dark Brown-Red
Protein by Biuret	40 - 100 mg/ml
WB-Cell Line/Tissue Extract	Rat Brain
Westen Blot-Titer	1:40,000-1:80,000
WB-Cell Line/Tissue Extract	NIH-3T3 Mouse Fibroblasts
WB-Westen Blot-Titer	1:20,000-1:40,000
Recommended Retest Period	-----
5 Years	

Specification: PRE.0.ZQ5.10000021007

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.